

(Deemed to be University Established Under Section 3 of UGC Act 1956) Pollachi Main Road, Eachanari Post, Coimbatore – 641 021. INDIA Phone: 0422-6471113-5, 6453777; Fax No: 0422-2980022-3 Email: info@karpagam.com; Web: www.kahedu.edu.in

### **DEPARTMENT OF BIOCHEMISTRY**

Subject	:	Practical VI- Clinical	Semester	:	III
		<b>Biochemistry and Animal</b>			
		Studies			
Subject code	:	17BCP312	Class	:	II M.Sc
					Biochemistry

### Lecture Plan

S. No	Duration of Period	Name of the Experiments	Support Materials
1	5	Estimation of glucose in serum	T1: 429-430
2	5	Estimation of cholesterol in serum	T1: 417-418
3	5	Estimation of urea in the urine and serum	T2: 124-125
4	5	Estimation of chloride in the urine and serum	T1: 454-455
5	5	Estimation of calcium in the urine and serum	T1: 50-51
6	5	Estimation of magnesium in the urine and serum	T1: 463-464
7	5	Analysis of urinary calculi	W1
8	5	Estimation of Bilirubin in serum (Kit method)	T2: 159-161

Prepared by Dr.P.Anusooriya, Department of Biochemistry, KAHE

Page 1/2

9	5	Estimation of triglyceride in serum (Kit method)	T2: 148-150
10	5	Estimation of HDL in serum (Kit method)	T1: 421-423
11	5	Handling of animals	4-6
12	5	Methods of injection	6-8
13	5	Induction of liver toxicity	W2
14	5	Assay of lipid peroxidation in rat liver	W3

### **Support Materials**

T1: Ramnick Sood, 2009. Medical Laboratory Technology (Methods and Interpretations), 6<sup>th</sup> Edition, Jaypee Brothers Medical Publishers (P) Ltd, New Delhi.

T2: Patel, A.H., 1994. Manual of Medical Laboratory Technology, 1<sup>st</sup> Edition, Navaneet Prakashan (L) Mumbai, India.

T3: Jann Hau, Steven J. Schapiro, 2010, Handbook of Laboratory Animal Science, 3<sup>rd</sup> Edition: Essential Principles and Practices, CRC Press.

W1: http://www.merckmanuals.com/professional/genitourinary-disorders/urinary-calculi/urinary-calculi

W2:http://www.hopkinsmedicine.org/healthlibrary/conditions/liver\_biliary\_and\_pancreatic\_disorders/alcoh ol-induced\_liver\_disease\_85,P00656

W3:http://www.researchgate.net/post/Protocol\_to\_analyze\_MDA\_TBARS\_method\_in\_homogenated \_samples



(Deemed to be University Established Under Section 3 of UGC Act 1956) Coimbatore - 641021. (For the candidates admitted from 2017 onwards) **DEPARTMENT OF BIOCHEMISTRY** 

SUBJECT	: Practical VI: C	linical Biochemistry and Animal Studies
SEMESTER	: III	
SUBJECT CODE	: 17BCP312	CLASS : II M.Sc.(BC)

#### **Course objectives**

- To know the analytical methods commonly used in the clinical laboratory.
- To know how can contribute the clinical laboratory to assess the health status of individuals.
- To understand the pathophysiology and molecular basis of the most prevalent diseases.

#### **Course outcome**

- Upon successful completion of this course, students will be able to:
- Explain the physiopathological bases and the biochemical markers of the most prevalent diseases in our population.
- Identify the principal analytical procedures used to measure biochemical magnitudes.
- Interpret and integrate the analytical data from the principal biochemical and molecular genetics tests for the screening, diagnosis, prognosis and monitoring of pathologies.

### **Clinical analysis**

- 1. Estimation of glucose in serum
- 2. Estimation of cholesterol in serum
- 3. Estimation of urea in the urine and serum
- 4. Estimation of uric acid in the urine and serum
- 5. Estimation of chloride in the urine and serum
- 6. Estimation of calcium in the urine and serum
- 7. Estimation of magnesium in the urine and serum
- 8. Analysis of urinary calculi
- 9. Estimation of Bilirubin in serum(Kit method)
- 10. Determination of A/G ratio
- 11. Estimation of triglyceride in serum (Kit method)
- 12. Estimation of HDL in serum (Kit method)

### **ANIMAL STUDIES (Group experiment)**

- 13. Handling of animals
- 14. Methods of injection
- 15. Induction of liver toxicity
- 16. Assay of lipid peroxidation in rat liver.

### REFERENCES

Jayaraman, J., (2007). Laboratory Manual in Biochemistry, New Age International Publishers New Delhi.

Sadasivam, S., and Manickam, A., (2009). Biochemical Methods, New Age International Publishers, New Delhi.

Singh, S.P., (2009). Practical Manual of Biochemistry, CBS Publishers, New Delhi.

Talib, V. H., (2003). A Handbook of Medical Laboratory Technology, CBS Publishers, New Delhi.



CLASS: II MSC BC COURSE NAME: Practical VI- Clinical Biochemistry & Animal studies COURSE CODE: 17BCP312 BATCH-2017-2019

### ESTIMATION OF GLUCOSE IN SERUM (O-Toluidine Method)

### AIM:

To estimate the amount of glucose present in the given serum and unknown solution.

### **PRINCIPLE:**

The aldehyde group of glucose condenses with o-toluidine in glacial acetic acid (aldehyde free colorless). On heating, giving an embral blue green N- glucosamine which is measured calorimetrically or photo metrically at 630nm. The intensity of color is proportion to the glucose concentration. The presence of thiourea stabilizes the o-toluidine reaction.



### **REAGENTS:**

✤ 10% TCA:

10 gm TCA dissolved in distilled water and made up to 100 ml.

**\*** O-toluidine reagent:

Dissolve 1.5gm pure thio urea in about 200 ml glacial aectic acid analytic grade, using gentle heat if necessary. Transfer to a 1 litre volumetric flask, washing in with glacial aetic acid. Add 60.0ml O-toluidine, mix, make to 1 liter with glacial acetic acid and mix by inversion. Keep in a brown bottle.

### Stock standard glucose solution:

200 mg of glucose is dissolved in 100 ml of distilled water.

### Working Standard (1mg/ml):

10 ml of stock standard was diluted to 100 ml of distilled water.

### **Procedure:**

In a serious of test tubes 0.5, 1.0, 1.5, 2.0, 2.5 ml of standard glucose solution corresponding the concentration 50, 100,150, 200, 250µg respectively.



### CLASS: II MSC BC COURSE NAME: Practical VI- Clinical Biochemistry & Animal studies COURSE CODE: 17BCP312 BATCH-2017-2019

- Serum Preparation: Pipette out 0.2 ml of the given serum sample and unknown solution were mixed with 1.8 ml of 10% TCA, Mix well, allow to stand 10 min at room temperature and centrifuge 10 minutes at 2500 rpm.
- > Transfer 1.0 ml of clear supernatant from step 2 in to the new test tube, labelled as test.
- ▶ In 2.5 ml of distilled water served as blank.
- > To all tubes add 5.0 ml of O-toluidine reagent, Mix by careful lateral shaking
- ▶ Heat in boiling water bath for 10 minutes.
- ➢ Cool in cold tap water for 4 min.
- ▶ Read OD of test and standard at 630 nm.

### Normal Value of Blood Glucose

Random	: 80-120 mg/dl
--------	----------------

Fasting : 60-100 mg/dl

Post prandial : Upto 120mg

### CALCULATION

### Unknown

The O.D value \_\_\_\_\_ corresponds to \_\_\_\_\_  $\mu g$  of glucose.

1.0 ml of unknown solution contains \_\_\_\_\_  $\mu$ g of glucose 100ml of unknown solution contains =  $\times$  100/1000

= \_\_\_\_ mg of glucose.

### Serum Sample

The O.D value \_\_\_\_\_ corresponds to \_\_\_\_\_  $\mu g$  of glucose.

1.0 ml of sample solution contains \_\_\_\_\_µg of glucose.

0.2 ml of sample solution contains =  $\frac{1.0 \times 0.2}{1.0 \times 0.2}$ 

100ml of sample solution contains =  $/1000 \times 100$ 

=\_\_\_\_ mg of glucose.



CLASS: II MSC BC COURSE NAME: Practical VI- Clinical Biochemistry & Animal studies COURSE CODE: 17BCP312 BATCH-2017-2019

	ESTIMATION OF GLUCOSE IN SERUM (O-toluidine Method)									
S.No	SOLUTION		Volume of	Volume of O-	Condition	Optical				
	Volume (ml)	concentration (µg)	distilled water (ml)	toluidine reagent (ml)		density 630nm				
Blank Standard	0.00	0.00	2.5							
S1	0.5	50	2.0							
S2	1.0	100	1.5							
S3	1.5	150	1.0		Heated in					
<b>S</b> 4	2.0	200	0.5	5ml	a boiling water bath					
S5	2.5	250	-		for 10					
Unknown					minutes					
U1	1.0	_	1.5							
U2	1.0	-	1.5							
Sample										
T1	1.0	-	1.5	↓						
T2	1.0	-	1.5							
	1	1	1	1						

### **Result:**

The amount of glucose present in the given serum and unknown solution is\_\_\_\_\_



CLASS: II MSC BC COURSE NAME: Practical VI- Clinical Biochemistry & Animal studies COURSE CODE: 17BCP312 BATCH-2017-2019

### ESTIMATION OF CHOLESTEROL BY ZAK'S METHOD

### AIM:

To estimate the amount of cholesterol in serum and unknown sample.

### **PRINCIPLE:**

Cholesterol in the presence of sulphuric acid it undergoes dehydration to form 3, 5cholestadiene. This is in turn oxidized and sulphonated to form red coloured cholestapolyene sulphonic acid in the presence of  $Fe^{3+}$  ions. The intensity of red color formed is proportional to the amount of cholesterol present in the serum. The colour intensity is measured by using a green filter 540 nm. This reaction has been employed by ZAK'S to estimate the cholesterol in serum and unknown sample.

### **REAGENTS:**

### **\*** Stock Standard Solution:

About 100 mg of cholesterol was dissolved and made up to 100 ml with glacial acetic acid.

### **\*** Working Standard:

About 10 ml of stock solution was made up to 100 ml with ferric chloride acetic acid reagent (concentration in 100  $\mu$ g/ml).

- ✤ Ferric chloride of 0.05% in glacial acetic acid
- \* Concentrated sulphuric acid
- ✤ Glacial acetic acid

### **PROCEDURE:**

- > 0.5 ml to 2.5 ml of working standard were pipetted out into clean test tubes.
- Serum Preparation: 0.1 ml of serum was taken and 9.9 ml of ferric chloride reagent was added. Mixed well and kept for 10 minutes at room temperature.
- > It was then centrifuged for 10 minutes at 3000 rpm.
- About of 0.5 ml and 1 ml of unknown and serum supernatant was taken in a test tubes. The volume was made up to 5.0 ml with ferric chloride and 3.0 ml of concentrated sulphuric acid were added.



### CLASS: II MSC BC COURSE NAME: Practical VI- Clinical Biochemistry & Animal studies COURSE CODE: 17BCP312 BATCH-2017-2019

- To prepare the blank, 5 ml of ferric chloride was mixed with 3 ml of concentrated sulphuric acid.
- > The test tubes were kept at room temperature for 15 minutes.
- The pinkish red colour formed was measured at 540 nm. Standard graph was drawn for the values obtained.
- From the standard graph the amount of cholesterol present in the unknown and serum sample can be calculated.

### Normal Value

Serum cholesterol: 150-200 µg/dl.

### CALCULATION

### Unknown

The O.D value \_\_\_\_\_ corresponds to \_\_\_\_\_ µg of cholesterol.

0.5 ml of unknown solution contains\_\_\_\_\_µg of cholesterol

100ml of unknown solution contains =  $\_$  × 100/0.5×1000

=\_\_\_\_ mg of cholesterol

### Serum Sample

The O.D value \_\_\_\_\_ corresponds to \_\_\_\_\_  $\mu g$  of cholesterol.

1.0 ml of sample solution contains  $\mu$ g of cholesterol.

10 ml of serum solution contains /0.5 x 10 µg of cholesterol

0.1 ml of serum solution contains\_\_\_\_\_µg of cholesterol

100ml of sample solution contains =  $(0.1 \times 100/1000)$ 

=\_\_\_\_ mg of cholesterol.

ESTIMATION OF CHOLESTEROL BY ZAK'S METHOD								
S.No	SOLUTION		Volume of	Volume of	Condition	Optical		
	Volume	concentration	FeCl <sub>3</sub> (ml)	Conc. H <sub>2</sub> SO <sub>4</sub>		density		

Exter | Ergipter | Exck Exter | Ergipter | Exck Example Constraints E

**KARPAGAM ACADEMY OF HIGHER EDUCATION** 

CLASS: II MSC BC COURSE NAME: Practical VI- Clinical Biochemistry & Animal studies COURSE CODE: 17BCP312 BATCH-2017-2019

	( <b>ml</b> )	(µg)		( <b>ml</b> )		540nm
Blank	0.00	0.00	5.0	♠		
Standard						
<b>S</b> 1	0.5	50	4.5			
S2	1.0	100	4.0			
<b>S</b> 3	1.5	150	3.5		Kept at	
<b>S</b> 4	2.0	200	3.0		room	
S5	2.5	250	2.5	3ml	temperature	
Unknown					for 15	
U1	0.5	_	4.5		minutes	
U2	0.5	_	4.5			
Sample						
T1	1.0	-	4.0			
T2	1.0	-	4.0	¥		
						1

### **Result:**

The amount of cholesterol present in the given serum and unknown solution is\_\_\_\_\_

# ESTIMATION OF UREA IN THE URINE AND SERUM (DAM- TSC METHOD) AIM:

To estimate the amount of urea present in blood and in urine.



CLASS: II MSC BC COURSE NAME: Practical VI- Clinical Biochemistry & Animal studies COURSE CODE: 17BCP312 BATCH-2017-2019

### **Principle:**

Urea reacts with diacetyl monoxime and thiosemicarbazide in the presence of ferric ions in a hot acidic medium to give a pink coloured compound. Ferric ions are provided by ferric chloride and acidic medium is provided by sulphuric acid and orthophosphoric acid present in acid reagent. The acid reagent is used to condense urea with diacetyl monoxime to yield a yellow coloured complex and further the colour is deepened to pink by the reaction with thioemicarbazide. The formed pink coloured product is spectrophotometrically measured at 540 nm.

### REAGENTS

### Stock Urea Standard Reagent:

Dissolve 1g pure dry urea and make up to 100 ml with distilled water. (1 ml= 10 mg)

### **Working Standard**:

Take 5 ml of stock solution and make up to 100 ml with distilled water  $(1ml = 500 \mu g)$ 

### **\*** DAM - TSC Reagent

Dissolve 100 mg of diacetyl monoxime, 20 mg of thiosemicarbazide, 900 mg Sodium Chloride in 100 ml of distilled water.

### Acid Reagent:

Add 1ml of orthophosphoric acid and 6 ml of concentrated sulphuric acid to about 75 ml of distilled water, cool the contents and add 0.1 ml of 10%  $FeCl_3$  solution. Dilute the solution to 100 ml with distilled water.

### **\*** Serum Preparation:

In to a dry test tube 3.4 ml of distilled water and 0.1 ml of blood. Mix well. Pipette 1.5 ml of 10% TCA and mix well. Centrifuged for 10 minutes.

### **PROCEDURE:**

- Pipette out 0.1 ml-0.5ml of the working standard solution (concentration varying from 50-250 µg) into a series of test tubes.
- Make up the given unknown sample to 100 ml with distilled water and pipette out 0.5 ml of this to a test tube.
- To all the test tubes, add distilled water to make up the volume to 0.5 ml. Take 0.5 ml of distilled water in another tube to serve as blank.



### CLASS: II MSC BC COURSE NAME: Practical VI- Clinical Biochemistry & Animal studies COURSE CODE: 17BCP312 BATCH-2017-2019

- Pipette out 3 ml of DAM-TSC Reagent and 3 ml of acid reagent into all the test tubes and keep at boiling water bath for about 15minutes and cool.
- Note the optical density of the pink colour developed colorimetrically at 540 nm. The intensity of the colour formed is directly proportional to the amount of urea present in it.
- Determine the concentration of the unknown sample from the graph and calculate the urea content in the given sample.

### Normal Value

Serum Urea: 10-50 mg/dL

Blood urea Nitrogen 5-25 mg/dL

### CALCULATION:

### Blood

Colorimetric reading \_\_\_\_\_ corresponds to \_\_\_\_\_ µg of urea.

0.5 ml of supernatant contains \_\_\_\_\_ µg of urea

100 ml of blood contains \_\_\_\_\_/1000 x100.

### Urine

Colorimetric reading \_\_\_\_\_ corresponds to \_\_\_\_\_ µg of urea.

1.5 ml of diluted urine contains \_\_\_\_\_ µg of urea

100 ml of diluted urine contains \_\_\_\_/1000 x100

1.0 ml of undiluted urine contain \_\_\_\_/1.0 x 1000 x 100

24hrs of urine sample contains\_\_\_\_/100 x 1500

ESTIMATION OF UREA DAM-TSC METHOD									
S.No	SOLUTION		Volum	Volum	Volum	Volum	Conditio	Optica	
	Volum	concentratio	e of	e of	e of	e of	n	1	
	e (ml)	n (µg)	Water	DAM	TSC	Acid		density	



CLASS: II MSC BC COURSE NAME: Practical VI- Clinical Biochemistry & Animal studies COURSE CODE: 17BCP312 BATCH-2017-2019

			( <b>ml</b> )	( <b>ml</b> )	( <b>ml</b> )	reagent		540nm
						(ml)		
Blank Standar	0.00	0.00	3.0	<b>f</b>	Ť	Ì		
d	0.1	50	2.5					
<b>S</b> 1	0.2	100	2.0					
<b>S</b> 2	0.3	150	1.5				Keep in	
<b>S</b> 3	0.4	200	1.0	1.0	1.0	3.0	boiling	
S4	0.5	250	0.5			5.0	water bath	
S5							for 30	
Sample	0.5	_	2.5				minutes	
T1	0.5	_	2.5					
T2								
Urine	1.5		1.5					
U1	1.5		1.5	+	+	↓		
U2								

### **RESULT:**

The amount of urea present in the given serum sample is  $\_\_\__mg/dL$ The amount of urea present in the given urine sample is  $\_\__g/24$  hrs urine.

### ESTIMATION OF CHLORIDE IN THE URINE AND SERUM (Van Slyke Method) AIM:

To estimate the amount of chloride present in urine and serum sample.

### **PRINCIPLE:**



### CLASS: II MSC BC COURSE NAME: Practical VI- Clinical Biochemistry & Animal studies COURSE CODE: 17BCP312 BATCH-2017-2019

This method uses the Volhard chloride estimation, protein being oxidized by the open Carius method using a solution of silver nitrate in concentrated nitric acid. The organic matter is destroyed and the chloride precipitated at the same the excess silver nitrate in then titrated with standard thiocyanate using ferric alum as the indicator.

### **REAGENTS**:

### > 0.05 N Silver nitrate

Dissolve 0.895g of silver nitrate in about 2 ml of water and then make up to 100 ml with concentrated nitric acid. The solution keeps in a brown bottle.

- > 0.02 N Thiyocyanate
- Concentrated nitric acid
- 5% Ferric alum

### PROCEDURE

Pipette a sample in to clean dry conical flask and run it slowly with constant shaking, 3ml of silver nitrate solution. Add 2 ml of concentrated nitric acid and heat over a Bunsen flame until digestion has occurred and the solution is pale Yellow in colour. If the colour is deep yellow it can be discharged by adding a few drops of saturated potassium permanganate solution and warming again. This should not be necessary with plasma. The digestion usually requires only one to two minutes. When the digestion is completed, cool and add 6 ml of 5% ferric alum. Then titrated with 0.02N thiocyanate until a reddish brown colour persisting for ten to fifteen seconds is obtained. To determine the standard add 2 ml of nitric acid and 6 ml of ferric alum solution to 3 ml of the silver nitrate solution, cool and titrate with 0.02N thiocyanate.

The difference between the two titration gives a measure of the amount of chloride in 1 ml serum and urine in terms of a 0.02 N solution.

CALCULATION: Serum Titration value of standard = \_\_\_\_ml Titre value of serum= \_\_\_ml



CLASS: II MSC BC COURSE NAME: Practical VI- Clinical Biochemistry & Animal studies COURSE CODE: 17BCP312 BATCH-2017-2019

The amount of chloride in 100 ml of serum = (\_\_\_ml titration of standard - \_\_\_ml titration of the unknown) x 20.

=\_\_\_\_mg of chloride

Urine

Titration value of standard = \_\_\_\_ml

Titre value of urine = \_\_\_\_ml

The amount of chloride in 100 ml of urine = (\_\_\_ml titration of standard - \_\_\_ml titration of the unknown) x 20.

=\_\_\_\_mg of chloride.

ESTIMATION OF CHLORIDE VAN SLYKE METHOD								
S. No	SOLUTION	Volume of Silver nitrate (ml)	Burette R Initial (ml)	Final (ml)	Volume of 0.02 N Potassium thiocyanate (ml)	Indicator		
1	Standard	3.0 3.0	0.0 0.0					
2	Serum 0.5 0.5	3.0 3.0	0.0 0.0			Ferric chloride		
3	Urine 1.0 1.0	3.0 3.0	0.0 0.0					

### **RESULTS:**

The amount of chloride present in 100 ml of serum is \_\_\_\_\_mg of chloride The amount of chloride present in 100 ml of urine is \_\_\_\_\_mg of chloride



CLASS: II MSC BC COURSE NAME: Practical VI- Clinical Biochemistry & Animal studies COURSE CODE: 17BCP312 BATCH-2017-2019

### ESTIMATION OF CALCIUM IN THE URINE AND SERUM

### AIM

To estimate the amount of calcium present in the given sample of urine and serum.

#### **PRINCIPLE:**

The calcium is precipitated as oxalate directly from the serum and after washing, the precipitate is dissolved in acid and titrated with permanganate.

### **REAGENTS:**

- 4% Ammonium oxalate
- > 2% ammonia
- > 0.01 N potassium permanganate

Prepare freshly before use by diluting stock 0.1 N solution

> Sulphuric acid

### **PROCEDURE:**

To 2 ml of serum add 2 ml of distilled water and 1 ml of 4% ammonium oxalate solution. Mix well and allow to stand overnight at room temperature. Tapered centrifuge tubes are more convenient for this determination. After precipitating the calcium, centrifuge and remove the supernatant fluid without disturbing the precipitate. With tapered tubes the supernatant fluid can be poured off. After doing so stand inverted on a filter paper and allow to drain for five minutes. Centrifuge again, pour off the supernatant fluid, repeat the draining on filter paper, wipe the mouth of the tube and add 2 ml of approximately normal sulphuric acid. Then placed the tubes in water bath for 70-80 °C for 5 minutes to dissolve the precipitate removed and titrated the contents hot with 0.01 N KMnO<sub>4</sub> solution in a burette till pale pink colour was got which persist for about a minute. Repeat the titration with a duplicates tubes.

Performed a blank titration with 2.0ml normal  $H_2SO_4$ . Kept in water bath for 5 minutes and titrated with KMnO<sub>4</sub> to a pink colour. The difference between the two tubes tire value gives the volume of 0.01 N KMnO<sub>4</sub> required to tire calcium oxalate.



CLASS: II MSC BC COURSE NAME: Practical VI- Clinical Biochemistry & Animal studies COURSE CODE: 17BCP312 BATCH-2017-2019

### CALCULATION:

### Urine

1 ml of 0.01 N KMnO<sub>4</sub> is equivalent to 0.2 mg of calcium. Hence since 2 ml of serum is used:

Mg. Calcium per 100 ml serum

= (Titration of unknown – titration of blank) x  $0.2 \times 100/2$ 

S. No	Volume of	Burette Reading		Volume of 0.02 N	Indicator
	Solution	Initial	Final	Potassium	
	(ml)	( <b>ml</b> )	( <b>ml</b> )	thiocyanate (ml)	
1	Blank				
	2.0	0.0			
	2.0	0.0			
2	Serum				$H_2SO_4$
	2.0	0.0			
	2.0	0.0			
3	Urine				
	2.0	0.0			
	2.0	0.0			

### **RESULT:**

The amount of calcium present in the 100 ml of urine is \_\_\_\_\_ mg of calcium.

# ESTIMATION OF MAGNESIUM IN THE URINE AND SERUM (Neill and Neely) AIM:

To determine the amount of Magnesium in given serum and urine sample.

### **PRINCIPLE:**

Exter ( Fragmer ) Brek Exter ( Fragmer ) Brek EXARPAGAM ACADEMY OF HIGHER EDUCATION (Deemed to be University) Prohobied Instruct and Table

# **KARPAGAM ACADEMY OF HIGHER EDUCATION**

### CLASS: II MSC BC COURSE NAME: Practical VI- Clinical Biochemistry & Animal studies COURSE CODE: 17BCP312 BATCH-2017-2019

Titan Yellow gives a red colour with magnesium Neill and Neely modified earlier methods by using gum ghatti as the colour stabilizer and including calcium in the standard. Calcium intensifies the colour so this allows for the effect of calcium present in the serum.

### REAGENTS

- Sodium tungstate, 10% solution
- Sulphuric acid 2/3N
- Gum gharrom 0.1 %: Suspend 0.1 g of powdered gum ghatti in a muslin bag in 100 ml of distilled water for twenty-four hours. The solution keeps well at room temperature.
- Titan Yellow, 0.05% solution: Dissolve 0.1 g of the powdered dye in 200ml of distilled water.

### Sodium hydroxide 4N

Stock standard solution containing 1mg /ml:

Dissolve 8.458 grams of MgCl<sub>2</sub>, 6H<sub>2</sub>O in distilled water and make up to a litre; or 10.094 grams of MgNH<sub>4</sub>PO<sub>4</sub>, 6H<sub>2</sub>O in 1N hydrochloric acid and make up to a litre, or dissolve 1.658 grams of reagent grade MgO heated to constant weight at red heat, in 10 to 20 ml. of concentrated hydrochloric acid and make up to a litre with water.

### Standard solution for use:

Dilute 1 ml of stock solution to 200 ml with water. This contains 5  $\mu$ g. Mg/ ml.

### > Calcium chloride solution, 0.05mg calcium/ml:

Dissolve 13.88 mf of Calcium chloride, (CaCl<sub>2</sub>) in water and make up to 100 ml.

### PROCEDURE

Dilute 1 ml of serum with 5 ml of distilled water and precipitate proteins by adding 2 ml of 10% sodium tungstate and 2 ml of 2/3 N sulphuric acid. Centrifuge to 5 ml of the supernatant fluid add in turn 1 ml distilled water, 1 ml of the gum ghatti, 1 ml of 0.05% Titan yellow and 2 ml of 4 N sodium hydroxide. At the same time put up 1 ml of calcium chloride and 5 ml of water and 1 ml of calcium chloride and 2.5 ml of the standard for use plus 2.5 ml of water as blank and standard respectively, completing these in the same way as the test. Read standard and unknown against the blank, using a green filter or with the instrument set at 520 millimicrons.

### CALCULATION



-

### CLASS: II MSC BC COURSE NAME: Practical VI- Clinical Biochemistry & Animal studies COURSE CODE: 17BCP312 BATCH-2017-2019

### Serum

Colorimeteric reading corresponds to µg of magnesium
0.2 ml of serum contains = $\_\ \mu g$ of magnesium
0.2 ml of sample was taken from 5 ml of the supernatant
The reaction mixture will contain
x 5/ 1000 x 0.2
1.0 ml serum contains mg of Magnesium
100 ml of serum contains x 100
mg of magnesium
Urine
Colorimeteric reading corresponds to µg of magnesium
1.0 ml of sample containing µg of magnesium
100 ml of urine sample contains = $x 100 / 1.0 \times 1000$
=mg of Magnesium

REAGENTS	Blank	<b>S1</b>	<b>S2</b>	<b>S</b> 3	<b>S4</b>	<b>S</b> 5	Ser	um	Ur	ine
							<b>S1</b>	<b>S2</b>	U1	U2
Volume of Working standard (ml)	-	0.2	0.4	0.6	0.8	1.0	0.2	0.2	0.2	0.2



CLASS: II MSC BC COURSE NAME: Practical VI- Clinical Biochemistry & Animal studies COURSE CODE: 17BCP312 BATCH-2017-2019

Concentration (µg)	-	1	2	3	4	5	-	-	-	-
Volume of water (ml)	2.5	2.3	2.1	1.9	1.7	1.5	2.3	2.3	2.3	2.3
Volume of CaCl <sub>2</sub> (ml)	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0
Volume of gum ghatti (ml)	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0
Volume of titan yellow (ml)	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0
Volume of NaOH (ml)	2.0	2.0	2.0	2.0	2.0	2.0	2.0	2.0	2.0	2.0
OD at 540 nm										

### RESULT

The amount of magnesium present in the given serum sample is \_\_\_\_\_ mg/ml.

### ANALYSIS OF URINARY CALCULI

### **PREPARATION OF REAGENTS**

- 1. 0.6 N HCl: 5 ml of conc. HCl in 100 ml of water
- 2. 1.7 M acetic acid: 10 ml of CH<sub>3</sub>COOH in 100 ml of water.
- **3. 0.02 M alcoholic p-Nitrobenzene azo resorcinol**: 0.5g in 100 ml ethanol.



CLASS: II MSC BC COURSE NAME: Practical VI- Clinical Biochemistry & Animal studies COURSE CODE: 17BCP312 BATCH-2017-2019

**4. Magnesium reagent:** 0.001% p-Nitrobenzene azo resorcinol in 2N NaOH stored in a polyethylene bolltle in a dark place.

5. Titan yellow: 7.5mg in 100 ml of polyvinyl alcohol

6. 0.3M lead acetate: 10g in 100 ml of water.

### **PROCEDURE:**

Washed the stones to free it from blood, mucous, etc., Stones submitted after sonic disintegration have blood and tissue adhering and or difficult to clean. Placed the stone in a beaker covered with several thickness of gauze held firmly in place with rubber band and washed under running water. Drained and removed the gauze carefully and dried the beaker and stones in an oven. Rinsed tiny crystal with water from squeeze bottle. Recorded the number, weight and dimensions of the stones. Describe briefly the colour and texture of exterior surface. The stones may be photographed for record purpose.

Cut saw or the stone was broken 50 as to examine the interior. Noted whether there was a foreign body which might have acted as a nucleus. A portion of each layer should be crushed and used for separate analysis to identify constituents.

Ashing the portion of the calculus in an open flame can be helpful in the identification of the constituents for flame analysis. To a small portion of the crushed stone was heated in a platinum dish until it glow. If the original bulk was not appreciably changed as ashing and if there was little or no darkening the calculus may be inorganic. If the same charged as was burned almost completely the calculus was found to be organic.

S. No	Experiment	Observation	Interface



CLASS: II MSC BC COURSE NAME: Practical VI- Clinical Biochemistry & Animal studies COURSE CODE: 17BCP312 BATCH-2017-2019

		Hard, smooth small or	Probably calcium oxalate
		Large small, round white or	Probably calcium
		Gray smooth, small, round	carbonate
	Dull stone smooth or rough	Probably uric acid	
1	1	White gray or yellow white, brown	Probably phosphate
1. Gross appearance	waxy appearance	Probably xanthine	
		Small, black, light weight	
		stone brittle when dry and	Probably fibrin
	Soft and plastic when wet waxy	Uro stealith	
		lustrous	Probably cysteine calculi
Chamic	al analysis of calculi		

### Chemical analysis of calculi

S No	Chamical group	<b>Descents and treatment</b>	Results and
5.110	Chemical group	Reagents and treatment	interpretation
		To a relatively large sample of	
1	Carbonata	pulverised stone in a test tube	Foaming effervescence
1.	Carbonate	added approximately 3 ml of 0.6	carbonate is present
		M HCl	
		Heated HCl solution from the	
		carbonate procedure in order to	
		solubilize the constituents. Cooled	
		and filtered via whatmann No. 1	If a white precipitate
2	Calcium oxalate	filter paper (5.5 cm dia) in to	forms calcium oxalate is
		another tube added 0.5 ml of	present
		saturated sodium acetate and	
		adjusted the pH to appr. 5 with 1.7	
		Molar CH <sub>3</sub> COOH	
Alterna	te pathway	·	·
1	For ovalate	Heated a small amount of	If effervescence occurs
1		pulverised stone in a porceleain	and there is no



			1
		dish. Cooled and placed in a test	effervescence in the
		tubes. Added approximately 2 ml	carbonate procedure
		of 0.6N HCl	oxalate is present.
			Moderate heat converts to
			carbonate with upon
			acidification yield CO <sub>2</sub>
			prolong intense heat will
			convert the oxalate in to
			oxide and invalidate this
			test
		To the acid extract added a pinch	
		of MnO <sub>2</sub> . The tubes should not be	Effervescence or tiny
		mixed or shaken. It may be	bubbles of $CO_2$ popping
2	Oxalate	necessary to warm the contents of	upward form the sediment
		the tubes very slightly to obtain the	of the tube is seen $MnO_2$
		reaction if only trace quantities of	oxidase the oxalate to
		oxalate are present.	$CO_2$ .
		If precipitate was formed with	If a white precipitate
		$CaC_2O_4$ procedure filtered using	forms no oxalate calcium
		Whatmann No. 1 filter paper.	is present.at pH of 5
		a) If no ppt is noted continued	oxalate combines with
		without filtering added	Ca <sub>2</sub> form a insoluble
2		0.5ml of 0.3 Molar	white precipitate of
3	Non oxalate calcium	potassium oxalate	$CaC_2O_4$
		(5.5g/100ml) adjusted the	
		pH to appr. 5	
		Another method	White precipitate or filim
		b) Acid extract on	reaction is most
		microscopic slide with 2 –	noticeable if the slide is
			1



-

	[				
		3 drops of saturated	heated over a dark		
		$(NH_4)_2C_2O_4$	background		
			Blue colour develops if		
		Dulyarized stope in a test tube and	present in trace amount		
		Furverised stone in a test tube and	this occurs upon standing		
4	Phosphates	one drop of 10% HCI and 2 drops	for a few minutes. The		
		of $H_2SO_4$ molibdate reagent and 1	blue colour is due to		
		drop of ANSA	reduced oxides of		
			molybdenum		
5			If the white precipitate		
		If precipitate was formed in the	forms magnesium and		
	Magnesium and phosphates	non-oxalate calcium produce	phosphate are present at		
		filtered through whatmann No.1	pH>8. The addition of		
		filter paper. if no precipitate was	NH <sub>4</sub> ions to solution		
		noted continued without filtrate	contain magnesium and		
		added NH <sub>4</sub> OH until a pH of 8 was	phosphate ions result in		
		obtained	the precipitate of NH <sub>4</sub>		
			MgPO <sub>4</sub>		
If no pr	ecipitate forms divided the	e solution into two portions and procee	eded as follows		
		To 1 portion added $0.5 \text{ ml of } 0.3 \text{M}$	If a white precipitate is		
6	Magnesium	Na <sub>2</sub> HPO <sub>4</sub> (5.5 $\alpha$ / 100ml)	formed magnesium is		
		1102111 04 (5.5g / 100111)	formed		
		To a $2^{nd}$ portion add 0.5 ml of 0.4	If a white precipitate is		
7	Phosphates	$M_{aso}$ (5 $\sigma/100$ ml)	formed phosphate is		
		Wi WigSO4 (3g/100mi)	present		
Alternate Procedure for Magnesium and phosphates:					
Heated a small amount of pulverized stone with 4 ml of 0.6N HCl. Cooled and divided the					
solution	n in to two parts. Analyzed	1 portion for magnesium and other for	or Phosphate.		
	a) Magnesium	Neutralized 1 portion with NH <sub>4</sub> OH			



		and added 1 ml of 0.02 M	Blue colour forms			
		alcoholic para nitro benzene azo	Magnesium is present			
		resorcinol (0.5g/100ml) ethanol				
		(or)				
		To the acid extract in the test tubes	Reddish purple reagent			
		added 2-3 drops of 20% NaOH	slowly becomes a definite			
		and 2-3 drops of magnesium	blue and (corn flour)			
		reagent	Precipitate forms.			
	b) Magnesium	Acid extract + 0.2 ml of Titan yellow (7.5mg/100ml) polyvinyl alcohol + 1 ml of 5N NaOH (2g/10ml)	If an orange red colour or precipitate forms magnesium is present			
8	Ammonium group	Acid extract was neutralized with 2-3 drops of 2.5 M NaOH + 0.5 ml Nesselr's reagent	Orange brown or rusty precipitate indicates ammonia positive result differentiate NH <sub>4</sub> urates from uric acid			
9	Cysteine a) Lead acetate Test	To an small amount of pulverized stone in a evaporating dish added 1-2 drops of 2.5M NaOH and heated. Added several drops of 0.3M lead acetate (10g/100ml) and heated again	Black precipitate is formed probably the present of cysteine is indicated			
Sulphic	les from other sources such	n as organic matter give a positive test	therefore. The presence of			
a cystei crystal.	a cysteine should be confirmed by the nitroprusside test or by a microscopic examination of crystal.					
	b) Nitroprusside	Boiled pulverized stone with 2 ml	Brick red colour was			
	test	of H <sub>2</sub> O. 2 ml of 1M NaCN	positive reaction which			



		(5g/100ml) waited for 5 minutes.	may fade to orange red
		Added 3 drops of a freshly	upon standing NaCN
		prepared solution of sodium nitro	converted cysteine to
		prusside (25 mg/ml) H <sub>2</sub> O	cysteine which then reacts
		(or)	with sodium forms a red
		Acid extract + 1 drop of 10%	colour
		NH <sub>4</sub> OH + 1 drop of 5% NaCN	
		waited for 5 minutes. Added 2-3	
		drops of sodium nitrioprusside	
		solution	
		Pulverized stones in a test tubes +	
		1 drop of 20% Na <sub>2</sub> CO <sub>3</sub> + 2 drops	Prompt deep blue colour
		of uric acid reagent	show uric acid
10	Urates and uric acids	(Or)	
		Acid extract was cooled and	
		neutralized with 2.5 M NaOH.	Pale blue colour is
		Added 1 ml of uric acid reagent	negative
Test for	l r differentiating between x	anthine and uric acid	
			The residue after
			evaporation is orange
		Pulverized stone in a test tube + 2	with uric acid but lemon
		drops of 10% HCl waited for 30	yellow with xanthine.
		minutes. Added 2 drops of 0.1 %	After addition of NaOH
11	Muroxide test uric acid	NaNO <sub>2</sub> waited for 30-60 minutes	an orange colour develops
	xanthine	added 2 drops of 0.5% NH <sub>4</sub>	becoming red with
		sulfamate + 2-3 drops sulpha dye	warming when xanthine is
		reagent	present a cherry red to
			purple colour develops
			immediately after addition



### CLASS: II MSC BC COURSE NAME: Practical VI- Clinical Biochemistry & Animal studies COURSE CODE: 17BCP312 BATCH-2017-2019

			of NaOH when uric acid
			is present
Bilary	calculi		
12	Cholesterol	To the pulverized stone in a test tube added 1 ml of CHCl <sub>3</sub> Heated over steam bath allowed to stand until the insoluble material settles, decanted into test tube. To this added 3 drops of acetic anhydride	If the test solution becomes red then blue and finally blue green in colour. It indicates the presence of cholesterol in
		and 1 drop of conc. H <sub>2</sub> SO <sub>4</sub> (Extraction can also be done with 3 ml of absolute alcohol)	the calculi
13	Bilirubin	Extracted a portion of the stone with several ml of ethanol added 0.5 ml of diazotised sulfanilic acid	Violet was appeared so bilirubin complexes with diazotized sulfanilic acid to form azobilirubin

### ANALYSIS OF URINARY CALCULI IN UNKNOWN SAMPLE

S.No.	Experiment	Observation	Inference
1.	To a relatively large sample of pulverized stone in a test tube added approximately 3.0ml of 0.6N HCl	Effervescence is observed	Presence of carbonate
2.	Heated the HCl solution from the carbonate procedure in order tom solubilize the constituents. Cooled & filteres through whatmann No:	White precipitated is formed	Presence of calcium oxalate.



	1 filter paper (5.5cm india) in			
	tube. Added 0.5 ml of saturated			
	sodium acetate and adjusted the			
	pH to approximately 5 with 1.7M			
	CH <sub>3</sub> COOH			
3.	a) Heated a small amount of			
	pulverized stone in a porceleain			
	dish, colled & placed the residue	Effervescence is		
	in a test tube added approximately	observed		
	2.0 ml of HCl.		Presence of oxalate	
	b) To acid extract added a pinch			
	of MnO <sub>2</sub> . The tube should not be	Tiny bubbles of CO <sub>2</sub>	Presence of oxalate is	
	mixed or shaken. It may be	popping upward from	confirmed	
	necessary to warm the content of	the sediment of the tube		
	the tube very slightly to obtain the	is observed		
	reaction if only trace quantities of			
	oxalate are present.			
	Pulverized stone in a test tube + 1			
4	drop of 10% HCl + 2 drops of	No blue colour was	Abaanaa of Dhaanhata	
4.	H <sub>2</sub> SO <sub>4</sub> molebdate reagent + 1	observed	Absence of Phosphate	
	drop of ANSA			
	Acid extract was neutralized with	Orange brown rusty		
5.	2-3 drops of 2.5M NaOH +0.5 ml	praginitate was formed	Presence of ammonia	
	of nessler's reagent	precipitate was formed		
	To a small amount of pulverized			
	stone in an evaporating dish	No black precipitate was		
6.	added 1 or 2 drops of 2.5M NaOH	formed	Absence of cysteine	
	and heated. Added several drops			
	of lead acetate (10mg/10g) and			



### CLASS: II MSC BC COURSE NAME: Practical VI- Clinical Biochemistry & Animal studies COURSE CODE: 17BCP312 BATCH-2017-2019

	heated again		
7.	Pulverized stone in a test tube + 1 drop of 20% Na <sub>2</sub> CO <sub>3</sub> + 2 drop of uric acid reagent	Prompt deep blue colour was observed	Presence of uric acid

### **RESULT:**

The given unknown salt contains carbonate, oxalate, ammonium and uric acid

ESTIMATION OF BILIRUBIN IN SERUM (DMSO/DIAZO Method)



CLASS: II MSC BC COURSE NAME: Practical VI- Clinical Biochemistry & Animal studies COURSE CODE: 17BCP312 BATCH-2017-2019

### AIM

To estimate the amount of bilirubin present in the given serum or plasma.

### PRINCIPLE

Sulfanillic acid reacts with sodium nitrite to form diazotized sulfanillic acid. Total bilirubin reacts with diazotized sulfanillic acid in the presence of DMSO to form azobilirubin, a red purple coloured product in acidic medium.

Bilirubin + Diazotized sulphanilic acid H<sup>+</sup> A zobilirubin H<sup>+</sup> Red purple colour complex

### SAMPLE

Serum (50µl is required)

Even though serum is recommended, plasma from blood treated with heparin can also be used for this test.

### **REAGENT COMPOSITION**

> Total bilirubin reagent

Sulfanilic acid

HCL

DMSO

Total bilirubin activator

Direct bilirubin reagent

Sulfanilic acid

HCL

- Direct bilirubin activator
- Bilirubin artificial standard

### PROCEDURE



### CLASS: II MSC BC COURSE NAME: Practical VI- Clinical Biochemistry & Animal studies COURSE CODE: 17BCP312 BATCH-2017-2019

S.No	Reagents	Total bilirubin		
		Sample Blank (µl)	Test (µl)	
1.	Total bilirubin reagent	1000	1000	
2.	Activator (T/D)	-	20	
3.	Serum	50	50	
Mix	well and incubate for exactly 5 minutes. I	Measure the absorbance of the	sample and standard	
	against the respective sa	mple blank at 546 or 532 nm.		

- > Take two test tubes Mark as sample blank and test and add 1000  $\mu$ l for both tubes.
- > Add 20  $\mu$ l activator for test alone.
- > Add 50  $\mu$ l of serum for all the tubes.
- Mix well and incubate for exactly 5 minutes.
- Measure the absorbance of the sample and standard against the respective sample blank at 546 or 532 nm.

### NORMAL RANGE

Total bilirubin : up to 1.2 mg/dl

Direct Bilirubin: up to 0.4 mg/dl

### **CALCULATION**

Total bilirubin = OD of test -OD of sample blank  $\times$  Factor

Total bilirubin concentration= OD of test –OD of sample blank

-----×10

OD of Standard

### RESULT

The amount of bilirubin present in the given serum sample is \_\_\_\_\_ mg/ml.

ESTIMATION OF TRIGLYCERIDS IN SERUM (GPO-PAP, End point Assay)



CLASS: II MSC BC COURSE NAME: Practical VI- Clinical Biochemistry & Animal studies COURSE CODE: 17BCP312 BATCH-2017-2019

### AIM

To estimate the amount of triglycerides present in the given serum or plasma.

### PRINCIPLE

Triglycerides are hydrolysed by lipoprotein lipase (LPL) to produce glycerol and fatty acids (FFA).In presence of glycerol kinase, ATP phosphorylates glycerol to produce glycerol-3-phosphate and ADP,glycerol-3-phosphate is further oxidized by glycerol-3-phosphate oxidase (GPO) to produce of peroxidase(POD)hydrogen peroxide couples with 4-amino anti pyridine (4 AAP) and 4-chlorophenol to produce red quinoneimine dye, absorbants of coloured dye is measured at 505 nm is proportional to TGL concentration in the sample.

### SAMPLE

Serum (10µl is required)

Serum or plasma is collected after 9-12 hours fasting.

### **REAGENT COMPOSITION**

- Triglyceride mono reagent
- > Triglyceride Standard

### PROCEDURE

Reagents	Blank (µl)	Standard (µl)	Test (µl)
Serum/Plasma	-	-	10
Reagent 2	-	10	-
Reagent 1	1000	1000	1000
Mix well and incubate at 37°C for 10 minutes. Measure the absorbance of the sample and standard			



### CLASS: II MSC BC COURSE NAME: Practical VI- Clinical Biochemistry & Animal studies COURSE CODE: 17BCP312 BATCH-2017-2019

against the respective sample blank at 505nm.

- > Take three test tubes mark as Standard, blank and test
- > Add 10  $\mu$ l of serum/plasma to test, 10  $\mu$ l of reagent 2 to standard.
- > 1000  $\mu$ l of reagent 1 to all the three tubes.
- > Mix well and incubate at  $37^{\circ}$ C for 10 minutes.
- Measure the absorbance of the sample and standard against the respective sample blank at 505 nm.

### NORMAL RANGE

Triglycerides concentration in mg/dl.

Normal: < 150 Border line high: 150-199 High: 200-499 Very high: > 500

### CALCULATION

Triglycerides (mg/dl)

Absorbance of test

-----×200

**RESULT** 

Absorbance of Standard

The amount of triglycerides present in the given serum sample is \_\_\_\_\_ mg/ml.

**ESTIMATION OF HDL IN SERUM (CHOD/POD METHOD)** 



CLASS: II MSC BC COURSE NAME: Practical VI- Clinical Biochemistry & Animal studies COURSE CODE: 17BCP312 BATCH-2017-2019

### AIM

To estimate the amount of HDL present in the given serum or plasma.

### PRINCIPLE

The reaction sequence employed in this assay is as follows:

Cholesterol ester + H<sub>2</sub>O  $\xrightarrow{CHE}$  Cholesterol + Free fatty acids Cholesterol + O<sub>2</sub>  $\xrightarrow{CHO}$  CHO Cholestenone + H<sub>2</sub>O<sub>2</sub> 2H<sub>2</sub>O<sub>2</sub> + Phenol+4- A mino antipyrine  $\xrightarrow{POD}$  Red Quinoneimine Complex + H<sub>2</sub>O

Cholesterol is determined after enzymatic hydrolysis and oxidation. Cholesterol esters are hydrolysed by the enzyme Cholesterol esterase to give free cholesterol and fatty acid molecules. This free cholesterol gets oxidized in the presence of Cholesterol oxidase to liberate Cholest-4 ene-3 one and peroxide. The indicator quinoneimine is formed from hydrogen peroxide and 4-aminoantipyrine in the presence of phenol and peroxidase. The intensity of this coloured complex is measured at 505 nm (500-540 nm) and is directly proportional to the cholesterol concentration present in the sample. On addition of the Precipitating reagent to the serum, followed by centrifugation, HDL fraction remains in the supernatant while the other Lipoproteins precipitate out.

### **SAMPLE PREPARATION**

Serum (10µl is required)

Take 0.5 ml of serum / plasma in to glass tube. Add 50  $\mu$ l precipitating reagent. Mix well, leave it at R.T for 10 minutes. Centrifuge at 3000 r.p.m for 10 minutes. Take the clear supernatant for HDL cholesterol estimation.

### **REAGENT COMPOSITION**

- Reagent 1:Cholesterol Enzyme reagent
- ▶ Reagent 2: Cholesterol Standard 200 mg/dl.
- Reagent 3:Cholesterol Precipitating reagent (Free)

### PROCEDURE



### CLASS: II MSC BC COURSE NAME: Practical VI- Clinical Biochemistry & Animal studies COURSE CODE: 17BCP312 BATCH-2017-2019

Reagents	Blank	Standard	Test
Enzyme reagent	1 ml	1 ml	1ml
Standard	_	10 µl	-
Supernatant sample	-	-	10 µl
Mix well and incubate at 37°C for 5 minutes. Measure the absorbance of the sample and standard			
against the respective sample blank at 510 nm.			

- > Take three test tubes mark as Standard, blank and test
- ➢ Add 1 ml reagent 1
- > Add 10  $\mu$ l of reagent 2 to standard and 10  $\mu$ l of supernatant.
- > Mix well and incubate at  $37^{\circ}$ C for 90 seconds.
- Measure the absorbance of the sample and standard against the respective sample blank at 510 nm.

### NORMAL RANGE

пл	cholesterol	
HDL	cnolesterol	

Men : 30-60 mg%

Women: 40-70 mg%

### CALCULATION

HDL cholesterol mg/dl =

Absorbance of HDL test

-----× 200

Absorbance of Standard

### RESULT

The amount of HDL cholesterol present in the given serum sample is \_\_\_\_\_ mg/ml.

### ANIMAL STUDIES



### CLASS: II MSC BC COURSE NAME: Practical VI- Clinical Biochemistry & Animal studies COURSE CODE: 17BCP312 BATCH-2017-2019

### HANDLING OF ANIMALS

Animal studies are really recognized fields in which animals are studied in a variety of cross displinary ways. Scholars from fields as diverse as art history, sociology, biology, physiology, literacy studies, geography, philosophy and ferminism, they seek to understand both human animal sections new and in the past and to understand animals as beings inn themselves separate from our knowledge of them.

#### Some experimental animals:

Mice, rat, Hamster, Guinea pig, Sheep, goat, Chicken, Dog, cat, most perfectly mices and rats are used.

#### **Rats:**

Rats should be acclimatized to handling to reduce stress. Rats should handle at the base of the tail using your fingers. Pick up the rats by placing the hand firmly over the back and the rib cage and restraining the head with thumb and fore finger immediately behind the mandibles. Presence of chassaignacs lumbercle, a thin plate of bone extending from the transverse process of the vertebrate makes it difficult to grasp the skin of the back. Holding the rat upside down keeps it distracted and reduces the chances of biting. A variety of resistant devices are available to assist handling the rats.

`Laboratory animals are inevitably subjected to human contact throughout their lives, during both husbandry and experiments. The use of appropriate and skilled handling is essential to ensure that animals readily accept or actively seek human contact and procedures are carried out efficiently. If routine handling procedures are aversive, animals are likely to develop anxiety and show exaggerated stress responses when approached. This is detrimental to animal welfare and will increase the difficulty of handling as animals attempt to avoid contact/restraint and may show defensive aggression. Handling stress can also be a major confounding variable and an unwanted source of variation within and between experiments. Good training in non-aversive handling has benefits for the animal, for the handler and for the reliability of data gained in experiments.

#### **Blood Sampling**



### CLASS: II MSC BC COURSE NAME: Practical VI- Clinical Biochemistry & Animal studies COURSE CODE: 17BCP312 BATCH-2017-2019

Collection of blood from laboratory animals is frequently necessary for a variety of experimental uses including determination of pharmacokinetics, antibody production, clinical pathology evaluation, etc. Blood may be collected from animals which are to survive the procedure or at sacrifice as a terminal event. Whereas there is no limitation on the amount of blood that may be collected terminally, the volume collected from animals surviving the collection is limited to prevent anemia and hypovolemia. As a general guide the 1-3-6 rule should be followed. The rule states that the average blood volume of most laboratory animals is 6% body weight (60 ml/kg); the most blood that can be reasonably expected from a terminal sacrifice is 3% body weight (30 ml/kg); and no more than 1% (10 ml/kg) body weight may be collected during any 2 week period from animals surviving the blood collection. Although venipuncture is generally a satisfactory method for survival blood collection of small quantities of blood.

### MICE

There are a variety of methods that are utilized to collect blood from mice. The techniques described below are recommended by OCV staff for survival (tail vein or orbital venous sinus) or terminal blood collection (cardiac).

### Lateral tail vein (Saphenous v.) venipuncture

The veins located on the lateral aspect of the mouse's tail are useful for collecting small volumes (< 0.1 ml) of blood. The technique for venipuncture is as described for IV injection except that a small volume of blood is aspirated into the syringe instead of injecting material. The use of a needle without a syringe, allowing the hub to fill with blood, and subsequent collection into a microhematocrit tube is useful when very small quantities of blood are needed.

### **Orbital venous sinus collection**

The sinus surrounding the globe of the mouse's eye is a useful site for collecting larger volumes of blood from surviving animals. The schematic provided in illustrates the location of the sinus. General anesthesia must be provided when collecting from this site. Under general anesthesia the mouse is grasped so that its back rests on the palm of your left hand (right hand if you are left-handed) with its head toward your thumb. The thumb is placed just just lateral to the animal's trachea so that the jugular vein on the same side as the eye from which you are collecting



### CLASS: II MSC BC COURSE NAME: Practical VI- Clinical Biochemistry & Animal studies COURSE CODE: 17BCP312 BATCH-2017-2019

blood is occluded and the fur on the animals head is drawn into the palm of your hand. This causes the animals eye to proptose (bulge) slightly. A 50 ul microhematocrit tube which has been broken in two is directed into the medial canthus of the eye rotating slightly as the tube is directed to a point directly behind the globe, inserting the non-broken end first. Sufficient pressure must be applied to cut through the fibrous layer which surrounds the sinus. Blood flows through the tube and occasionally around the tube once the sinus has been penetrated. After blood collection, the tube is removed and the eyelids closed and a dry cotton pledget is applied over the eye with gentle pressure to prevent retroorbital hemorrhage. In general blood should not be collected from the same eye more than 3 times, allowing at least 1 week between collections. An antibiotic opthalmic ointment may be applied following bleeding.

### Cardiac puncture (diaphragmatic approach)

Cardiac puncture is the preferred technique for terminal collection of large blood volumes. General anesthesia is administered and the animal placed on a solid surface with its ventrum exposed The xyphoid process is palpated at the caudal aspect of the animal's sternum. A notch is present on both sides of this process. Negative pressure should be applied, by placing slight backward pull on the plunger, once it has been inserted beneath the skin. Reflux of blood is apparent once the needle has penetrated the heart.

### RATS

### Lateral tail vein (Saphaneous v.) vein puncture

The procedure for collecting blood from the rat's tail vein is similar to the technique described for the mouse. A slightly smaller gauge needle (24 ga or larger gauge) can be utilized. Because of the vein's size, larger blood volumes (approximately 1 ml) may be obtained from adult rats.

### **Orbital venous plexus**

The technique describing blood collection from the mouse's orbital venous sinus should be followed for orbital venous plexus collection in the rat. The only difference in that the vessels surrounding the rat's globe are a network of small veins rather than a blood filled sinus and the fibrous connective tissue surrounding the plexus is quite dense. Therefore, the broken end of the hematocrit tube which serves as a cutting edge should be inserted into the plexus. Remember this



CLASS: II MSC BC COURSE NAME: Practical VI- Clinical Biochemistry & Animal studies COURSE CODE: 17BCP312 BATCH-2017-2019

technique must be performed under general anesthesia and post bleeding hemostasis is essential to prevent complications.

### Cardiac puncture (diaphragmatic approach)

The technique for cardiac puncture from the rat is identical to that described for the mouse except a longer and smaller gauge needle is recommended. A 5 - 10 cc syringe should be used if large blood volumes are desired. This procedure is performed as a terminal event only and general anesthesia is required. The animal must be sacrificed at the completion of the procedure prior to awakening from anesthesia.

### **METHODS OF INJECTION**

### RATS

#### Subcutaneous

SC injections are performed in rats using the same technique as was described for mice with the following differences. The volume of material administered can be increased to approximately 5 ml per site in an adult rat (>300 grams). Syringe size should be increased proportionately and needles should be 22 ga or larger gauge. May be performed in any area of loose skin along the back or flank. Tenting the skin between the shoulder blades or over the rump creates an appropriate pocket for injection. Inserting the needle under the skin along the flank results in the outline of the needle clearly visible when correctly situated.



#### Intramuscular



### CLASS: II MSC BC COURSE NAME: Practical VI- Clinical Biochemistry & Animal studies COURSE CODE: 17BCP312 BATCH-2017-2019

IM injections may be performed in the rat. Injection volumes are limited to 0.25 ml site because of limited muscle mass. Either the quadriceps muscles located on the cranial aspect of the femur or the hamstrings on the caudal aspect of the femur can be used. Care must be taken to avoid depositing material on or near the sciatic nerve which runs along the caudal aspect of the femur in the thigh. The needle is directed through the skin into the muscle belly approximately 3-4 mm. Aspiration should be attempted before injecting to determine that accidental penetration of a blood vessel has not occured.



### Intravenous

IV injection technique for the rat is similar to the mouse. However, the vessels are more difficult to visualize, especially in adult rats. The skin overlying the vessels in adults becomes quite thick, making vascular access much more difficult. For this reason the preferred site for vascular access is near the tail base. Injection volumes administered to an adult rat should not exceed 2 mls and large volumes should be administered slowly to avoid vascular overload. The technique describing IV administration and needle size in mice should be followed. Dilate the blood vessels by warming the mouse. Once warmed, place the mouse in a restrainer. Insert the needle as low as possible towards the tip of the tail, since the vein is very superficial at the tip. The vein will clear from the injection site to the base of the tail if properly situated, whereas ballooning around the injection site will occur if the needle is not properly seated you will note ballooning at the injection site.



CLASS: II MSC BC COURSE NAME: Practical VI- Clinical Biochemistry & Animal studies COURSE CODE: 17BCP312 BATCH-2017-2019



### Intraperitoneal

The technique for IP injections in rats is virtually identical to mice. Rats should be restrained with their abdomen exposed and their head held downward. The injection site, method and needle size is as described for mice. Because of their larger size < 5.0 mls of material can be administered to an adult rat. Holding the mouse in dorsal recumbency, insert the needle in a position below the bend of the knees; left or right of the midline. Avoid the midline to prevent penetrating the bladder. Angle the needle approximately  $45^{\circ}$  to the body.



### ASSAY OF LIPID PEROXIDATION IN ETHANOL INTRODUCED LIVER TOXICITY IN RATS

### Introduction of liver toxicity

Rats are given 40% ethanol (v/v 2.0 ml/100 g body weight) for 30 days.



CLASS: II MSC BC COURSE NAME: Practical VI- Clinical Biochemistry & Animal studies COURSE CODE: 17BCP312 BATCH-2017-2019

### Principle

In this method, Malondialdehyde and other thiobarbituric acid reactive substance (TBARS) were measured by their reactivity with thiobarbituric acid in the acidic condition to generate a pink colour chromophore which was measured at 535nm.

#### Reagents

- ✤ 0.67 % Thio barbituric acid
- ✤ 10% TCA
- Tris HCL buffer

#### Procedure

The tissue homogenate was prepared in Tris HCL buffer. 1 ml of liver homogenate as combined with 2ml of TBA-TCA-HCL reagent and mixed thoroughly. The mixture was kept in boiling water bath for 15 minutes. After cooling the tubes they were centrifuged at 100 g for 10 minutes and supernatant was taken for measurement. A series of standard solution were treated in a similar mannar. The absorbance of chromophore was measured at 535 nm against a reagent blank that contain no tissue homogenate.

TBARS were expressed as n moles of MDA formed/ mg protein.