

19MBU511A

Semester – V  
(4H – 2C)

## MANAGEMENT OF HUMAN MICROBIAL DISEASES - PRACTICAL

Instruction Hours / week: L: 0 T: 0 P: 4

Marks: Internal: 40 External: 60 Total: 100

End Semester Exam: 9 Hours

### COURSE OBJECTIVES

- To provide a strong base in the fundamentals of pathogens.
- To learn techniques and methods used in the cultivation and isolation of pathogens.
- To obtain with the knowledge about the habitat and characteristics of pathogens.

### COURSE OUTCOME

Involves the identification, classification, and characterization of pathogenic species.

### EXPERIMENTS

1. Diagnosis of respiratory tract disease.
2. Diagnosis of urinary tract disease.
3. Diagnosis of gastrointestinal tract disease.
4. Identification of dermatophytes

### SUGGESTED READINGS

1. Ananthanarayan R. and Paniker C.K.J. (2009) Textbook of Microbiology. 8<sup>th</sup> edition, University Press Publication
2. Brooks G.F., Carroll K.C., Butel J.S., Morse S.A. and Mietzner, T.A. (2013) Jawetz, Melnick and Adelberg's Medical Microbiology. 26<sup>th</sup> edition. McGraw Hill Publication.
3. Goering R., Dockrell H., Zuckerman M. and Wakelin D. (2007) Mims' Medical Microbiology. 4<sup>th</sup> edition. Elsevier.
4. Willey JM, Sherwood LM, and Woolverton CJ. (2013) Prescott, Harley and Klein's Microbiology. 9<sup>th</sup> edition. McGraw Hill Higher Education.
5. Madigan MT, Martinko JM, Dunlap PV and Clark DP. (2014). Brock Biology of Microorganisms. 14<sup>th</sup> edition. Pearson International Edition.

## **Experiment No: 1**

### **DIAGNOSIS OF RESPIRATORY TRACT DISEASE**

**INTRODUCTION:** The respiratory tract is far more prone to infection than other organ systems, as it frequently comes into contact with numerous pathogens present throughout the environment. Infections of the respiratory tract are caused by a wide range of viruses and bacteria that may vary seasonally and can impact the upper and lower respiratory tract, affecting different patient populations with varying degrees of severity.

**AIM:** To diagnose the received specimen and to isolate and to identify the organism suspected for respiratory tract disease.

**SAMPLE:** Sputum is the most usual specimen. Apart from it bronchial washings, biopsies, gastric aspirate, CSF, pleural fluids were also collected in certain cases.

Sputum is a mixture of bronchial secretion and inflammatory exudates coughed out and expectorated in some infections. During bacterial infection sputum will be purulent and contain yellow, green and opaque materials. Usually disposable wide mouthed plastic container of about 100 ml capacity is used to collect specimen.

### **SPECIMEN PROCESSING**

**Homogenization:** The specimen is usually processed in 2 ways:

1. By means of adding sterile glass beads.
2. Equal volume of 1% buffered pancreatin was added and shaken well and was placed in water bath for an hour. The suspension was shaken well every 15 minutes and this helps in the uniform distribution of the specimen.

### **PROCEDURE**

#### **Microscopy**

**Gram staining:** Make the smear and perform the gram staining to see the morphology.

**Acid fast staining (or) Ziehl Nelson staining:** The heat fixed smear of specimen was flooded with carbolfuchsin and heated. Then the smear was washed with water and decolorized with acid alcohol (20% H<sub>2</sub>SO<sub>4</sub> and then with 95% ethanol) for 2 minutes. The smear was then counterstained with Loeffler's methylene blue and then observed under oil immersion.

#### **Culture method**

**Lowenstein Jensen's medium:** Inoculate the specimen onto two bottles of LJ medium and incubated at 37°C and inspected for growth of typical colonies weekly once.

**Blood agar:** Prepare Blood agar plates and streak the given specimen on the surface of the medium and incubate at 37 °C for 24 h.

**Nutrient agar:** Prepare the nutrient agar plates and streak the given specimen on the surface and incubate at 37 °C for 24 h.

**Macconkey agar:** Prepare Macconkey agar plates and streak the given specimen on the surface and incubate at 37 °C for 24 h.

**Mannitol salt agar:** Prepare Mannitol salt agar plates and streak the given specimen on the surface and incubate at 37 °C for 24 h.

### **Biochemical characteristics**

#### **Catalase test**

Media-Nutrient agar

Culture – 24 h old culture

Reagent - H<sub>2</sub>O<sub>2</sub> (3%)

#### **Coagulase test**

**Tube coagulase:** It is used to detect free coagulase. Add 0.4 ml of (1.5 dilution) plasma to 0.1 ml of overnight broth and incubate for 3-6 h and observe the tubes for coagulation.

**Slide coagulase:** It is used to detect bound coagulase. Place 2 drops of saline onto slide. Take a loopful of inoculum and introduce and following which add 2 drops of undiluted plasma. Then observe for coagulation.

**Gelatin liquefaction test:** Prepare Nutrient Gelatin media and inoculate with the given specimen. Incubate the plates at 37°C for 24 h.

### **INTERPRETATION OF RESULTS**

#### **Microscopy**

##### **Gram staining:**

##### **Acid fast staining (or) Ziehl Nelson staining:**

#### **Culture method**

##### **Lowenstein Jenson's medium:**

##### **Blood agar:**

##### **Nutrient agar:**

##### **Macconkey agar:**

##### **Mannitol salt agar:**

#### **Biochemical characteristics**

##### **Catalase test:**

**Coagulase test:**

**Tube coagulase:**

**Slide coagulase:**

**Gelatin liquefaction test:**

Based on the observation of the microscopy, cultural and biochemical characteristics, the possible organism in the specimen may be interpreted as \_\_\_\_\_.

## **Experiment No. 2**

### **DIAGNOSIS OF URINARY TRACT DISEASE**

**INTRODUCTION:** Urinary tract infections occur more frequently in women than men due to shortness of female urethra. *E.coli* is the principle urinary tract pathogen. Other less common cause of UTI is *Staphylococci*, *Enterococci*, *Proteus*, *Klebsiella* and other *Coliforms* and *Pseudomonas aeruginosa*.

Most infections originate from the patients bowel flora (endogenous infection). Infections are usually caused by the entry of microorganisms from external sources. Microorganisms at the opening of the urethra, careless personal hygiene and sexual intercourse contribute to the incidence of UTI. Opportunistic organisms present in intestine can also cause the infection.

**AIM:** To diagnose the specimen and to isolate and to identify the organism suspected for urinary tract disease.

**SPECIMEN COLLECTION:** Specimen should be collected carefully without any urethral contamination.

**Clean catch midstream urine:** The first portion of the urine which is more likely to be contaminated must be avoided as the bacteria will have multiplied to high levels after overnight incubation in the bladder. Clean catch of midstream urine was collected in sterile wide mouth screw cap tubes.

## **PROCEDURE**

### **Microscopic examination**

**Wet mount method:** Place a drop of fresh uncentrifuged urine on slide and cover with cover glass and examine under light microscope.

**Gram staining and Motility test:** Place Gram stained smear of urine on slide and cover with cover glass and examine under light microscope. Carry out the Hanging drop technique using the specimen to check the motility.

**Culture methods:** Dilute the urine specimen and streak on nutrient agar, macconkey agar, mannitol salt agar, cystine–lactose–electrolyte-deficient agar (CLED) plates and incubate the agar plates at 37 °C for 24 h.

### **Biochemical characteristics**

#### **Indole test**

Media- peptone broth

Reagent – 0.2 ml of Kovac's reagent

Incubation time & temperature – 37 °C for 24 h

### **Methyl red test**

Media- MRVP broth

Reagent – methyl red

Incubation time & temperature – 37 °C for 48 h.

### **Voges-Proskauer test**

Media- MRVP broth

Reagent – Barritt's reagent

Incubation time & temperature – 37 °C for 48 h.

### **Citrate test**

Media- Simmon's citrate agar

Reagent – bromothymol blue

Incubation time & temperature – 37 °C for 48 h.

### **Carbohydrate fermentation test**

Media- Peptone broth

Reagent - bromothymol blue

Sugars employed – 1 % (w/v) Glucose, Sucrose, Lactose, Mannitol

Incubation time & temperature – 37 °C for 24-48 h.

### **Urease test**

Media-Nutrient agar

Reagent - Christensen's Urea Agar

Incubation time & temperature – 37 °C for 24 h.

### **Catalase test**

Media-Nutrient agar

Reagent - H<sub>2</sub>O<sub>2</sub> (3%)

Culture – 24 h old culture

## **INTERPRETATION OF RESULTS**

### **Microscopic examination**

#### **Wet mount method:**

#### **Gram staining and Motility test:**

#### **Culture methods:**

#### **Biochemical characteristics**

#### **Indole test**

#### **Methyl red test**

#### **Voges-Proskauer test**

## **Citrate test**

## **Carbohydrate fermentation test**

## **Urease test**

## **Catalase test**

Based on the observation of the microscopy, cultural and biochemical characteristics, the possible organism in the specimen may be interpreted as \_\_\_\_\_.

### **Experiment No. 3**

#### **DIAGNOSIS OF GASTROINTESTINAL TRACT DISEASE**

**INTRODUCTION:** Gastrointestinal infections are viral, bacterial or parasitic infections that cause gastroenteritis, an inflammation of the gastrointestinal tract involving both the stomach and the small intestine. Symptoms include diarrhea, vomiting, and abdominal pain. Diarrhea is a condition where intestinal discharge is expelled in a liquid state which is known as the stool. Diarrhea is accompanied by discomfort, loss of water and electrolytes and it leads to severe dehydration. Pathogenic organisms such as *Vibrio cholera* causes to severe dehydration, *E. coli* causes travelers' diarrhea, some intestinal viruses and Rota viruses can also cause diarrhea.

**AIM:** To diagnose the specimen and to isolate and to identify the organism suspected for gastrointestinal tract disease.

**SPECIMEN:** Stool sample

**COLLECTION AND TRANSPORT:** The specimen was collected in a clean, dry, disinfectant free wide mouthed container without urine contamination. When fecal sample is not available rectal swab is preferable. It can transport using sterile Cary Blair medium.

#### **PROCEDURE**

##### **Microscopic examination**

**Gram staining & Motility test:** Make the smear and perform the gram staining to see the morphology. Carry out the Hanging drop technique using the specimen to check the motility.

**Culture methods:** Inoculate the specimen on various media like EMB agar, Macconkey agar, thiosulphate citrate bile salt media, xylose lysine deoxycholate agar, SS agar and incubate for 16-24 h at 37 °C.

##### **Biochemical tests**

###### **Indole test**

Media- peptone broth

Reagent – 0.2 ml of Kovac's reagent

Incubation time & temperature – 37 °C for 24 h

###### **Methyl red test**

Media- MRVP broth

Reagent – methyl red

Incubation time & temperature – 37 °C for 48 h.

###### **Voges-Proskauer test**

Media- MRVP broth

Reagent – Barritt's reagent



Incubation time & temperature – 37 °C for 48 h.

**Citrate test**

Media- Simmon's citrate agar

Reagent – bromothymol blue

Incubation time & temperature – 37 °C for 48 h.

**Catalase test**

Media-Nutrient agar

Reagent - H<sub>2</sub>O<sub>2</sub> (3 %)

Culture – 24 h old culture

**INTERPRETATION OF RESULTS**

**Microscopic examination**

**Gram staining & Motility test:**

**Culture methods**

**EMB agar:**

**MacConkey agar:**

**Thiosulphate citrate bile salt media:**

**Xylose lysine deoxycholate agar:**

**SS agar:**

**Biochemical tests**

**Indole test**

**Methyl red test**

**Voges-Proskauer test**

**Citrate test**

**Catalase test**

Based on the observation of the microscopy, cultural and biochemical characteristics, the possible organism in the specimen may be interpreted as \_\_\_\_\_.

## **Experiment No. 4**

### **IDENTIFICATION OF DERMATOPHYTES**

**INTRODUCTION:** Dermatophytes are a unique group of closely related filamentous fungi that invade keratinized cutaneous structures, including the stratum corneum, nails, and hair, of humans and animals, resulting in an infection referred to as dermatophytosis, ringworm, or tinea. Diagnosis of dermatomycoses is based on the detection of septate hyphae by direct microscopic examination of clinical samples. Direct microscopic examination is rapid and inexpensive but does not provide genus or species identification. *Arthroderma* now contains 21 species, *Ctenomyces* one species, *Epidermophyton* one species, *Lophophyton* one species, *Microsporum* three species, *Nannizzia* nine species and *Trichophyton* 16 species. In addition, two new genera have been introduced: *Guarromyces* containing one species and *Paraphyton* three species. Although the number of genera has increased, the species that are relevant to routine diagnostics now belong to smaller groups, which should enhance their identification.

**AIM:** To identify the dermatophytes

**SPECIMEN:** Skin, hair and nail samples

**COLLECTION AND TRANSPORT:** Skin, hair and nail samples were collected from clinically suspected cases of dermatophytosis under aseptic conditions using 70 % alcohol into a sterile black paper.

### **PROCEDURE**

**KOH wet mount:** The specimens collected were subjected to KOH wet mount preparation. 10 % KOH used for skin and hair samples and kept for 10-15 min and 40 % KOH used for nail samples, kept overnight. A small amount of sample is taken and added to the drop of KOH placed on a glass slide and covered with a glass slip and observed under microscope.

**Culture:** The specimen was inoculated onto Sabouraud's dextrose agar with 0.05 % Chloramphenicol and 0.5 % Cycloheximide and incubated at 28 °C for up to four weeks, and was observed periodically for growth and (reverse) pigmentation.

**Tease mount by lactophenol cotton blue:** A small portion of a colony was picked and suspended in two drops of lactophenol cotton blue placed on a clean slide. The mycelial mat was teased apart with dissecting needles, covered with cover-slip and observed under microscope for presence of aseptate slender hyphae, macro and microconidia and their arrangement.

**Dermatophyte Test Media (DTM):** All isolated dermatophytes were inoculated onto DTM and incubated at 28 °C for 7 days and observed for color change.

**Urease test:** All isolated dermatophytes were inoculated onto Christensen's Urea Agar and incubated at 28 °C for 7 days and observed for colour change. This test is to differentiate between *T. mentagrophytes* and *T. rubrum*.

## INTERPRETATION OF RESULTS

**KOH wet mount:**

**Culture:**

**Tease mount by lactophenol cotton blue:**

**Dermatophyte Test Media (DTM):**

**Urease test:**

Based on the observation of the microscopy, cultural and biochemical characteristics, the possible organism in the specimen may be interpreted as \_\_\_\_\_.